


EU-SolidAct
Centralisation of Analyses - Visit Baseline (D1), D3, D8 and D15
(Naso-oro-pharyngeal samples)

Analyses	Sampling Materials	Transfer Materials	Storage Conditions	Recipient	Transport Conditions	Carrier
Naso-oro-pharyngeal swab	1 sampling flocked swab + Transport medium 3ml	2x 2ml cryotubes	-80°C	Theradis Pharma	With Dry ice	Specialised carrier



Sampling Kit



Naso-oro-pharyngeal swab + Transport medium 3ml



Cryotubes

- Use the sampling kit identified: « **Visit D3 / Naso-oro-pharyngeal swab** » (Example)
- Verify the expiry date noted on the sampling kit - **Never use a kit which has expired !**
- Prepare and identify all the necessary materials prior to each visit

- Complete the labels for the sampling tubes and the cryotubes



Theradis Pharma



Study File

- **Take the patient requisition form** from the sampling kit corresponding to the visit
- **Complete** the sections reserved for the investigator and sampler
- Send the Patient Requisition Form (white copy) **by email or by fax** to Theradis Pharma **within 24 hours**:

E-mail: solidact@theradispharma.com
Fax: + 33 (0)4 97 10 08 78
- **Keep the pink copy** of the patient requisition form in the study file

